

# Importing and Processing data with Python

## 1 Introduction

This guide outlines the essential Python libraries and steps for importing and processing EEG data, especially for applications involving ERP (Event-Related Potential) analysis. We will cover data import, filtering, and segmentation, using libraries such as MNE, SciPy, NumPy, and Matplotlib.

## 2 Required Libraries for EEG Data Analysis

The following libraries are commonly used in EEG data analysis:

- **MNE (Machines for Neural Experiments)**: A Python package specifically for EEG data analysis, offering functions for preprocessing, filtering, artifact removal, time-frequency analysis, and visualization.
- **SciPy**: Provides scientific routines and enables loading of MATLAB files for EEG data.
- **NumPy**: Supports multi-dimensional arrays, matrices, and mathematical functions, essential for data manipulation.
- **Matplotlib**: A visualization library used to plot EEG signals, spectrograms, and more.

## 3 Importing EEG Data

When loading EEG data from a `.mat` file, the contents are stored in a dictionary-like object with keys representing different data segments. EEG data is usually organized as a multi-dimensional array with the `shape` attribute denoting the size of each dimension:

```
Number of channels = eeg.data.shape[0]  
Number of time samples = eeg.data.shape[1]  
Number of trials = eeg.data.shape[2]
```

Using the MNE library, we create an `Info` object to store details like channel names, sampling frequency, and channel types. The sampling rate should be set according to the data, e.g., 512 Hz for BCI Competition III data.

## 4 Filtering EEG Data

EEG signals span a range of frequencies (typically 1-30 Hz for relevant activity), but also contain noise:

- **Low-frequency noise:** Caused by slow movements or wire shifts.
- **High-frequency noise:** Caused by muscle movements or electromagnetic interference.

To retain only the relevant band (1-30 Hz), we apply a band-pass filter:

- **Low-pass filter:** Removes high-frequency noise.
- **High-pass filter:** Removes slow drifts.

**Tip:** Always filter the data before segmenting into shorter epochs. Use MNE's `.filter()` method as follows:

```
raw_filt = raw.copy().filter(low_cut, high_cut)
```

where `low_cut` is typically set to 0.1 Hz and `high_cut` to 30 Hz.

## 5 Segmentation into ERP Epochs

After filtering, segment the EEG data into **epochs**, which are time-locked to specific events (e.g., stimuli or responses). The `events` array and an event dictionary (mapping codes to labels) are used to define these segments:

```
events, events_dict = mne.events_from_annotations(raw_filt)
```

The `events` array has three columns:

- **First column:** Index of each event in the data.
- **Second column:** Usually zero, ignored.
- **Third column:** Event code as an integer.

### 5.1 Event Dictionaries

MNE assigns unique integer values to each event code and stores this mapping in the event dictionary. This helps standardize event codes from different data sources.

## 6 References

- Neurotist. (2024, February 12). *Importing MATLAB Files into Python: A Step-by-Step Guide for EEG Data Analysis with MNE*. Medium. Retrieved from <https://medium.com/@neurotist/importing-matlab-files-into-python-a-step-by-step-guide>
- *Filtering EEG Data - Neural Data Science in Python*. (n.d.). Retrieved from [https://neuraldatascience.io/7-eeg/erp\\_filtering.html](https://neuraldatascience.io/7-eeg/erp_filtering.html)